


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Resorcinol - ReagentPlus[®], 99%

Product Number: 307521
CAS Number: 108-46-3
MDL: MFCD00002269
Formula: C₆H₆O₂
Formula Weight: 110.11 g/mol



TEST

Specification

Appearance (Color)	White to Pink
Appearance (Form)	Conforms to Requirements
Powder, Chips, Granules or Crystals	
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.5 %
Solubility (Turbidity)	Clear
50mg/ml (5%), H ₂ O	
Solubility (Color)	Colorless to Very Faint Yellow
Recommended Retest Period	-----
9 Years	

Specification: PRD.0.ZQ5.10000019768

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.