

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Hexyl methacrylate - 98%, contains 100 ppm MEHQ

Product Number:

462373

CAS Number:

142-09-6

MDL:


MFCD00015283

Formula:

C₁₀H₁₈O₂

Formula Weight:

170.25 g/mol



TEST

Specification

Appearance (Color)

Colorless to Faint Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 97.5 %

Stabilizer

Approximately 100ppm MEHQ

Specification: PRD.1.ZQ5.10000005287

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

