

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Human IgG1-Biotin antibody produced in mouse - clone 8c/6-39, purified immunoglobulin, buffered aqueous solution

Product Number:

B6775

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Form)

Liquid

Appearance (Turbidity)

Clear to Slightly Hazy

Appearance (Color)

Colorless

IgG Concentration (OD 280nm)

2.0 - 4.0 mg/ml

ELISA - Titer

1:2000-1:4000

Recommended Retest Period

3 Years

Specification: PRE.0.ZQ5.10000017956

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.