

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum(IV) chloride - $\geq 99.99\%$ trace metals basis**Product Number:****379840****PtCl₄**

CAS Number:

13454-96-1

MDL:

MFCD00011182

Formula:

Cl₄Pt

Formula Weight:

336.89 g/mol

TEST**Specification**

Appearance (Color)

Faint Orange to Very Dark Brown

Appearance (Form)

Conforms to Requirements

Powder and/or Chunks

ICP Major Analysis

Confirmed

Confirms Pt Component

Trace Metal Analysis

 ≤ 100.0 ppm

Purity

Meets Requirements

 $\geq 99.99\%$ Based On Trace Metals Analysis

Specification Date : 02/14/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.