

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Formamide - ReagentPlus[®], ≥99.0% (GC)

Product Number:

F7503

CAS Number:

75-12-7

MDL:


MFCD00007941

Formula:

CH₃NO

Formula Weight:

45.04 g/mol



TEST

Specification

Appearance (Turbidity)

Clear

Appearance (Color)

Colorless to Very Faint Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 99.0 %

Specification: PRD.3.ZQ5.10000004320

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.