

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
SSPE Buffer 20× Concentrate – for Northern and Southern blotting, solution

Product Number: S2015

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
pH	7.2 - 7.6
Neat at room temperature	
DNAse, Exonuclease Detection	None Detected
NICKase,Endonuclease Detection	None Detected
RNAse Detection	None Detected
Suitability	Suitable
Suitable for Southern Transfer. Six or seven bands detected.	

Specification: PRD.0.ZQ5.10000064921

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.