


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Maleic acid disodium salt hydrate - ≥98%

Product Number: M9009
CAS Number: 25880-69-7
MDL: MFCD00150133
Formula: C₄H₂Na₂O₄ · xH₂O
Formula Weight: 160.04 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
50 mg/mL, H ₂ O	
Water (by Karl Fischer)	≤ 19 %
Sodium (Na) anhydrous	26 - 30 %
Titration with H ₂ SO ₄	≥ 98 %
Anhydrous	

Specification: PRD.0.ZQ5.10000030377

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.