

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification


Product Name:  
Oxalyl chloride solution - 2.0 M in methylene chloride

Product Number: **310670**

MDL: MFCD00000704

Formula: C<sub>2</sub>Cl<sub>2</sub>O<sub>2</sub>

Formula Weight: 126.93 g/mol



### TEST

### Specification

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Purity Titration by AGNO<sub>3</sub>

1.95 - 2.10 M

Specification: PRD.0.ZQ5.10000020121

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.