

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Silver tetrafluoroborate - 98%

Product Number: 208361
CAS Number: 14104-20-2
MDL: MFCD00003408
Formula: AgBF₄
Formula Weight: 194.67 g/mol

AgBF₄

TEST	Specification
------	---------------

Appearance (Color)	White to Off-White
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Titration by KSCN	54.0 - 56.8 %
% Ag	

Specification: PRD.0.ZQ5.10000025580

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.