

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium hexafluorophosphate - ≥99%

Product Number:

200913

CAS Number:

17084-13-8

MDL:

MFCD00011412

Formula:

F6KP

Formula Weight:

184.06 g/mol



TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Powder or Crystals or Chunk(s) or

Powder with Chunk(s) or Crystals with Chunk(s)

ICP Major Analysis

Confirmed

Confirms K and P Components

Purity

Confirmed

> = 99 %

Specification: PRD.3.ZQ5.10000004707

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.