

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4-(Trifluoromethyl)phenyl isocyanate - 99%

Product Number:

250821

CAS Number:

1548-13-6

MDL:

MFCD00002028

Formula:


C₈H₄F₃NO

Formula Weight:

187.12 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 98.5 %

Specification: PRD.0.ZQ5.10000026403

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.