

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:
Pyrrolidine - ≥99%, FGProduct Number: **W352316**

CAS Number: 123-75-1

Formula: C₄H₉N

Formula Weight: 71.12 g/mol



TEST	Specification
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.441 - 1.445
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 99.0 %
Specific Gravity at 25 deg C	0.856 - 0.860
Color Test	≤ 300 APHA
Arsenic (As)	≤ 3.0 ppm
Cadmium (Cd)	≤ 1.0 ppm
Mercury (Hg)	≤ 1.0 ppm
Lead (Pb)	≤ 10.0 ppm
Expiration Date Period 5 Years	-----

Specification: PRD.5.ZQ5.10000009966

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.