

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Tetraamminepalladium(II) nitrate solution - 10 wt. % in H₂O, 99.99%

Product Number:

423823



MDL:

MFCD00011591

Formula:

H₁₂N₆O₆Pd

Formula Weight:

298.55 g/mol

TEST	Specification
Appearance (Color)	Yellow
Appearance (Form)	Liquid
Gravimetric Analysis % Pd; With Dimethylglyoxime	2.8 - 4.3 %
ICP Major Analysis Confirms Palladium Component	Confirmed
Purity 99.99% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 150.0 ppm

Specification: PRD.0.ZQ5.10000039673

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.