

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)


Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
21-Deoxycortisol

Product Number: **P9521**  
CAS Number: 641-77-0  
Formula: C<sub>21</sub>H<sub>30</sub>O<sub>4</sub>  
Formula Weight: 346.5 AMU



| TEST                           | Specification                 |
|--------------------------------|-------------------------------|
| Appearance (Form)              | Powder                        |
| Appearance (Colour)            | White to White w/ Yellow Cast |
| Solubility (Solvent)           | Ethanol                       |
| Solubility (Conc)              | 9.80 - 10.20 mg/ml            |
| Solubility (Heating)           | Yes                           |
| Solubility (Sonication)        | Yes or No                     |
| Solubility (Turbidity)         | Clear                         |
| Solubility (Color)             | Colorless                     |
| Elemental Anal. (%C anhydrous) | 72.31 - 73.31 %               |
| Purity (TLC)                   | ≥ 98.00 %                     |
| NMR (Solvent)                  | DMSO-d <sub>6</sub>           |
| Identity by NMR                | Consistent                    |
| Recommended Retest Period      | -----                         |
| 3 Years                        |                               |

Specification: PRE.0.ZQ5.10000002285

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.