

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Acetone oxime – Vetec™ reagent grade, 98%

Product Number:

V900596

CAS Number:


127-06-0

Formula:

C₃H₇NO

Formula Weight:

73,09 g/mol



TEST

Specification

Appearance (Colour)

White to Yellow

Appearance (Form)

Crystals

Infrared spectrum

Conforms to Structure

Solubility (Turbidity)

Clear to Slightly Hazy

c = 100 mg/ml (10%); Water

Solubility (Colour)

Colorless to Yellow

Purity (GC)

≥ 97.5 %

Specification: PRE.0.ZQ5.10000021081

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.