

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

**Product Number:** 635065  
**CAS Number:** 1317-80-2  
**MDL:** MFCD00011269  
**Formula:** O<sub>2</sub>Ti  
**Formula Weight:** 79.87 g/mol



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Crystals or Crystals with Chunk(s) or Crystalline Chunks
Particle Size	Confirmed
10mm x 10mm x 0.5 mm	
Orientation 110	Confirmed
ICP Major Analysis	Confirmed
Confirms Titanium Component	
Purity	Meets Requirements
> = 99.9% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 1000.0 ppm

Specification: PRD.1.ZQ5.10000032829

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.