

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Carbon nanotube, single-walled - >95% (carbon as SWCNT)

Product Number: **775533**
CAS Number: 308068-56-6



TEST	Specification
------	---------------

Appearance (Color)	Black
Appearance (Form)	Powder
Purity	≥ 95
Carbon Content	
Sheet Resistance	≤ 600
Ohms/sq at 85% VLT	
Registered Trademark	Confirmed
SWeNT is a Registered Trademark of Chasm Advanced Materials	

Specification: PRD.2.ZQ5.10000031835

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.