


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(R_p)-1-[(1S)-(1-Aminoethyl)]-2-(diphenylphosphino)ferrocene - 97%

Product Number: **697354**
CAS Number: 607389-84-4
Formula: C₂₄H₂₄FeNP
Formula Weight: 413.27 g/mol



TEST

Specification

Appearance (Color)	Yellow to Orange to Brown
Appearance (Form)	Conforms to Requirements
Crystals or Powder	
Infrared spectrum	Conforms to Structure
Carbon	67.4 - 72.2 %
Optical Rotation	322.0 - 393.0 deg
c = 0.1%, CHCl ₃	
Proton NMR spectrum	Conforms to Structure
Supplier information	Conforms
Product of Kanata Chemical Technologies, Inc.	

Specification: PRD.1.ZQ5.10000061594

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.