

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium Acetate Buffer Solution – 100 mM, pH 5.0, 0.2 µM filtered

Product Number:

SAE0157

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)

Colorless to Almost Colorless

Appearance (Form)

Liquid

Appearance (Turbidity)

Clear

pH

4.9 - 5.1

Specification: PRD.0.ZQ5.10000102636

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

