


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Carvacrol - ≥98%, FCC, FG

Product Number: **W224502**
 CAS Number: 499-75-2
 MDL: MFCD00002236
 Formula: C₁₀H₁₄O
 Formula Weight: 150.22 g/mol



TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.521 - 1.526
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.0 %
Residual Solvent	≤ 100 ppm
Toluene	
Specific Gravity	0.974 - 0.980
Solubility (Turbidity)	Clear to Slightly Hazy
1ml/4ml 60% Ethanol	
Solubility (Color)	Colorless to Faint Yellow
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.3.ZQ5.10000031424

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

