

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification


Product Name:
L-Carvone - natural, 99%, FG

Product Number: **W224908**

CAS Number: 6485-40-1

Formula: C₁₀H₁₄O

Formula Weight: 150.22 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.495 - 1.502
Infrared spectrum	Complete
Angular Rotation	-62.0 - -57.0 °
Purity (GC)	≥ 98.5 %
Arsenic (As)	≤ 3 ppm
Mercury (Hg)	≤ 1 ppm
Cadmium (Cd)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Kosher Status	Confirmed
Expiration Date Period	-----
60 Months	

Specification: PRD.0.ZQ5.10000027321

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.