

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

tert-Butyl 1-(Boc-amino)-3-cyclopentene-1-carboxylate - 97%

Product Number:

683000

CAS Number:

521964-59-0

MDL:


MFCD08460299

Formula:

C₁₅H₂₅NO₄

Formula Weight:

283.36 g/mol

**TEST****Specification**

Appearance (Color)

White to Tan

Appearance (Form)

Crystal to Powder

Infrared spectrum

Conforms to Structure

Carbon

61.4 - 65.8 %

Nitrogen

4.7 - 5.1 %

Specification: PRD.0.ZQ5.10000054929

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.