

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
trans-Ferulic acid - 99%

Product Number: **128708**
CAS Number: 537-98-4

Formula: C₁₀H₁₀O₄
Formula Weight: 194.18 g/mol



TEST	Specification
------	---------------

Appearance (Colour)	Off White to Yellow and Faint Beige to Beige
Appearance (Form)	Powder or Powder with Chunk(s)
Titration with NaOH	98.5 - 101.5 %
Purity (GC)	≥ 98.5 %
Infrared Spectrum	Conforms to Structure

Specification: PRE.0.ZQ5.10000017254

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.