

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

5-Chlorothiophene-2-sulfonyl chloride - 96%

Product Number:

544272

CAS Number:

2766-74-7

MDL:

MFCD00051667

Formula:


C₄H₂Cl₂O₂S₂

Formula Weight:

217.09 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Appearance (Color)

Colorless to Yellow to Brown

Appearance (Form)

Solid or Liquid

Proton NMR Spectrum

Conforms to Structure

Purity (GC)

≥ 95.5 %

Specification: PRD.2.ZQ5.10000018415

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.