

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

**Product Number:** P8291  
**CAS Number:** 7757-79-1  
**MDL:** MFCD00011409  
**Formula:** KNO<sub>3</sub>  
**Formula Weight:** 101.10 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Solubility (Turbidity)	Clear
20g/150ml in Water	
Solubility (Color)	Colorless
Titration	≥ 99.0 %
Purity by Titration (Dry Basis)	
Loss on Drying	≤ 1.00 %
Plant Cell Culture Test	Pass
Cell Culture Test	Pass
Recommended Retest Period	-----
3 Years	

Specification: PRD.3.ZQ5.10000033894

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.