

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

(S)-4-[(1-Methylheptyloxy)carbonyl]phenyl 4'-octyloxy-4-biphenylcarboxylate - 99%

Product Number:

408859

CAS Number:

112901-67-4

MDL:


MFCD00192300

Formula:

C36H46O5

Formula Weight:

558.75 g/mol



TEST

Specification

Appearance (Color)

White

Appearance (Form)

Powder or Chunks

Infrared spectrum

Conforms to Structure

Carbon

76.2 - 78.5 %

Optical Rotation

21.0 - 27.0 deg

c = 1%; Chloroform

Purity (HPLC)

≥ 98.5 %

Specification: PRD.0.ZQ5.10000011845

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.