


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
8-Hydroxyquinoline – crystalline

Product Number: H6878
CAS Number: 148-24-3
Formula: C₉H₇NO
Formula Weight: 145.16 AMU



TEST	Specification
Appearance (Color)- White to light brown	Conforms
Appearance (Form)- Crystals or Powder	Conforms
Purity	98.5 - 101.5 %
Identification IR	Conforms
Expiration Date Period 2 Years	-----

Specification: PRE.1.ZQ5.1000020115

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.