


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Chloro(1,5-cyclooctadiene)rhodium(I) dimer - 98%

Product Number: 227951
CAS Number: 12092-47-6
MDL: MFCD00012415
Formula: C₁₆H₂₄Cl₂Rh₂
Formula Weight: 493.08 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Orange, Yellow-Orange or Orange-Brown	
Appearance (Form)	Conforms to Requirements
Powder, Crystals or Crystalline Powder	
Infrared spectrum	Conforms to Structure
Carbon	38.0 - 39.9 %

Specification: PRD.0.ZQ5.10000053877

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.