


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium thiosulfate - ReagentPlus®, 99%

Product Number: 217263
CAS Number: 7772-98-7
MDL: MFCD00003499
Formula: Na₂O₃S₂
Formula Weight: 158.11 g/mol



| TEST | Specification |
|------|---------------|
|------|---------------|

| | |
|---------------------------------------|-------------------------------|
| Appearance (Color) | White |
| Appearance (Form) | Powder or Crystal or Chunk(s) |
| Titration by Iodine | 98.5 - 101.5 % |
| ICP Major Analysis | Confirmed |
| Confirms Sodium and Sulfur Components | |
| Solubility (Turbidity) | Clear |
| 100 mg/ml (10%), H ₂ O | |
| Solubility (Color) | Colorless |
| Recommended Retest Period | ----- |
| 3 Years | |

Specification: PRD.0.ZQ5.10000023967

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.