


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
N,N-Diglycidyl-4-glycidyoxyaniline

Product Number: 412813
CAS Number: 5026-74-4
MDL: MFCD00192036
Formula: C₁₅H₁₉NO₄
Formula Weight: 277.32 g/mol
Storage Temperature: 2 - 8 °C



| TEST | Specification |
|------|---------------|
|------|---------------|

| | |
|----------------------------|------------------------------|
| Appearance (Color) | Colorless to Yellow to Brown |
| Appearance (Form) | Liquid |
| Infrared spectrum | Conforms to Structure |
| Purity (GC) | ≥ 90.0 % |
| Epoxide Equivalent (g/mol) | 95 - 106 |
| Viscosity | 550 - 950 cps |
| At 25 Degrees Celsius | |

Specification: PRD.0.ZQ5.10000046881

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.