

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Calcium chloride - anhydrous, beads, - 10 mesh, ≥99.9% trace metals basis

Product Number: 449709 CaCl₂
CAS Number: 10043-52-4
MDL: MFCD00010903
Formula: CaCl₂
Formula Weight: 110.98 g/mol

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Beads
ICP Major Analysis Confirms Ca Component	Confirmed
Trace Metal Analysis	≤ 1000.0 ppm
Purity ≥99.9% Based On Trace Metals Analysis	Conforms

Specification: PRD.1.ZQ5.10000027786

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.