


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
cis-3-Hexenyl hexanoate - FG

Product Number: **W340308**
 CAS Number: 31501-11-8
 MDL: MFCD00036552
 Formula: C₁₂H₂₂O₂
 Formula Weight: 198.3 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Refractive index at 20 °C	1.435 - 1.439
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.0 %
Color Test	≤ 100 APHA
HPLC (% Minor 1) (alpha-tocopherol)	≤ 0.1 %
Acid Value	≤ 2.0 ml
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Organoleptics	Pass
Expiration Date Period 5 Years	-----

Specification: PRD.2.ZQ5.10000005639

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

