


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

2,2'-Azobis(2-methylpropionitrile) solution 0.2 M in toluene

Product Number:

714887

Formula:

C₈H₁₂N₄

Formula Weight:

164.21 g/mol

TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Proton NMR spectrum	Conforms to Structure
Proton NMR	0.10- 0.30 M

Remarks:

Specification Date : 08/16/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.