

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-UCP-1 antibody produced in rabbit – affinity isolated antibody, buffered aqueous solution

Product Number:

U6382

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	1.0 - 1.5 mg/ml
WB-Cell Line/Tissue Extract	Rat Brown Fat Mitochondria
WB Titer	1:1000-1:2000
IHC - Tissue	Mouse Brown Fat
IHC - Fixation	FFPE+ Antigen Retrieval by Decloaking
IHC - Titer	1:500-1:1000
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000016971

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.