

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Fomblin[®] Y - LVAC 14/6, average mol wt 2,500

Product Number: **317934**
CAS Number: 69991-67-9



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Viscosity	120.00 - 160.00
Kinematic Viscosity, Centistokes	
mg KOH/g - Neutralization Numb	≤ 0.03
%O (Mass)	
Evap. Weight Loss (22H/149°C)	≤ 5.00 %
Vapor Pressure (100 °C)	≤ 2.70 Pa
Vapor Pressure (20 °C)	≤ 2.70 Pa

Specification: PRD.0.ZQ5.10000030842

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.