


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Fluorometholone - ≥98%

Product Number: **F9381**
CAS Number: 426-13-1
MDL: MFCD00056461
Formula: C₂₂H₂₉FO₄
Formula Weight: 376.46 g/mol



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
25 mg/mL, Pyridine	
EmM	14.8 - 15.6
Lambda Max 239nm, MeOH	
Purity (HPLC)	≥ 98 %

Specification: PRD.2.ZQ5.10000040328

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.