

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 716928

Formula Weight: 196.97 g/mol

Storage Temperature: 2 - 8 °C

Au

TEST	Specification
Appearance (Color)	Grey to Black
Appearance (Form) Colloid	Conforms
Size Axial Diameter	22.5 - 27.5 nm
Size Longitudinal Size	65.7 - 80.3 nm
Concentration Concentration of Gold	45.0 - 55.0 ppm
Miscellaneous Assay Concentration of Palladium	35.0 - 45.0 ppm
Recommended Retest Period 1 Year	-----

Specification: PRD.2.ZQ5.10000036048

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.