


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium hydrogencarbonate - -40-+140 mesh, ≥95%

Product Number: 401676
CAS Number: 144-55-8
MDL: MFCD00003528
Formula: CHNaO3
Formula Weight: 84.01 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Titration by HCL	≥ 95.0 %
Impurity	2.0 - 5.0 %
% Na ₂ CO ₃	
X-Ray Diffraction	Conforms to Structure
Size	Conforms
Typically ca. 90% within -40 to + 140 mesh	

Specification: PRD.0.ZQ5.10000033493

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.