

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Sodium methoxide solution - 25 wt. % in methanol

Product Number:

156256



MDL:

MFCD00012179

Formula:

CH₃NaO

Formula Weight:

54.02 g/mol

TEST

Specification

Appearance (Color)

Colorless to Very Faint Green-Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Titration with HCl

23.5 - 26.5 %

Specification: PRD.2.ZQ5.10000032401

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.