

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
β-Glucuronidase from Escherichia coli - Type IX-A, lyophilized powder, 1,000,000-5,000,000 units/g protein
(30 min assay)

Product Number: G7396
CAS Number: 9001-45-0
MDL: MFCD00130629

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Appearance (Color)	White to Light Yellow to Light Brown
Appearance (Form)	Powder
% Protein (Biuret)	30 - 70
units/g Protein	1000000 - 5000000
One Sigma or Modified Fishman unit will liberate 1.0 microgram of phenolphthalein from phenolphthalein glucuronide per hour at pH 6.8 at 37 deg C (30 minute assay).	
Recommended Retest Period	-----
2 years	

Specification: PRD.1.ZQ5.10000029596

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.