

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Angiotensin Converting Enzyme from porcine kidney - lyophilized powder, ≥ 10 units/mg protein (Bradford)

Product Number: **A2580**
CAS Number: 9015-82-1
MDL: MFCD00130325

Storage Temperature: -20 °C

TEST	Specification
------	---------------

units/mg protein

 ≥ 10

Unit definition: One unit will produce 1.0 umole of hippuric acid from Hippuryl-His-Leu per min in 50 mM HEPES and 300 mM NaCl at pH 8.3 at 37 °C

Specification: PRD.0.ZQ5.10000030517

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.