


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-(Diethylamino)propylamine - $\geq 99\%$

Product Number: **549975**
CAS Number: 104-78-9
MDL: MFCD00008218
Formula: C₇H₁₈N₂
Formula Weight: 130.23 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared Spectrum	Conforms to Structure
Water (by Karl Fischer)	$\leq 0.5 \%$
Color Test (APHA)	≤ 10
Hazen	
Assay	$\geq 99.0 \%$
Purity	
Note	Confirmed
Product of Arkema Inc.	

Specification: PRD.1.ZQ5.1000009167

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.