

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium 2-phenoxyphenyltrifluoroborate - 97%

Product Number:

766615

CAS Number:

213745-39-2

Formula:

C₁₂H₉BF₃KO

Formula Weight:

276.10 g/mol

TEST

Specification

Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder or Chunks
Fluorine NMR Spectra	Conforms to Structure
Boron HNMR Spectra	Conforms to Structure
Proton NMR spectra	Conforms to Structure
Carbon	50.4 - 54.0 %
Purity Titration NAOH by HCL	96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000025265

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.