

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
17-Dimethylaminoethylamino-17-demethoxygeldanamycin

Product Number: D5193
CAS Number: 467214-20-6
Formula: C₃₂H₄₈N₄O₈
Formula Weight: 616.75 g/mol
Storage Temperature: -20 °C

TEST	Specification
------	---------------

Appearance (Color)	Purple
Appearance (Form)	Powder
Solubility (Color)	Purple
Solubility (Turbidity)	Clear
10 mg/ml, EtOH	
Purity (TLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000027759

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.