


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
D-2-Aminobutyric acid - 98%

Product Number: 116122
CAS Number: 2623-91-8
MDL: MFCD00064414
Formula: C₄H₉NO₂
Formula Weight: 103.12 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form) Crystalline Flakes	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Assay	≥ 97.5 %
Optical Rotation c = 5%, 6N HCl or 2%, 5N HCl	-21.0 - -19.0 deg

Specification: PRD.2.ZQ5.10000046835

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.