

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Fluorophenylhydrazine hydrochloride - 97%

Product Number:

153974

CAS Number:

2924-16-5

MDL:


MFCD00012934

Formula:

C₆H₇FN₂ · HCl

Formula Weight:

162.59 g/mol



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Conforms to Requirements

Fibers, Powder or Crystals

Infrared spectrum

Conforms to Structure

Titration by AgNO₃

96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000042045

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.