

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
TiterMax® Classic Adjuvant - liquid

Product Number: **H4397**

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Emulsification

Pass

Forms emulsion with 0.85% NaCl

Note

Conforms

TiterMax Classic Adjuvant contains three essential ingredients: a block copolymer, CRL89-41, squalene (a metabolizable oil) and a unique microparticulate stabilizer.

TiterMax is a registered trademark of TiterMax USA, Inc.

Specification: PRD.0.ZQ5.10000033731

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.