

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Potassium isopropenyltrifluoroborate - 97%

Product Number:

723916

CAS Number:

395083-14-4



Formula:

C3H5BF3K

Formula Weight:

147.98 g/mol

TEST

Specification

Appearance (Color)	White to Beige
Appearance (Form)	Powder or Solid
Proton NMR spectrum	Conforms to Structure
Fluorine NMR Spectra	Conforms to Structure
Boron HNMR Spectra	Conforms to Structure
Carbon	23.5 - 25.2 %
Purity Titration NAOH by HCL	96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000021748

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.