

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Polypropylene Isotactic, average M_w ~190,000, average M_n ~50,000

Product Number:

427896


CAS Number:

9003-07-0

MDL:

MFCD00084447

Formula:

C₃H₆**TEST****Specification**

Appearance (Color)

White or Colorless

Appearance (Form)

Conforms to Requirements

Pellets or Translucent Beads

Infrared spectrum

Conforms to Structure

Melt Index: g/10 minutes

28.0- 40.0

230 Degrees Celsius, 2.16kg

Remarks:

Specification Date : 11/29/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.