


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium hexahydroxyplatinate(IV)

Product Number: 464708
CAS Number: 12325-31-4
MDL: MFCD00014236
Formula: H6Na2O6Pt
Formula Weight: 343.1 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Yellow to Orange
Appearance (Form)	Powder
ICP Major Analysis	Confirmed
Confirms Sodium and Platinum Components	
Gravimetric Analysis	51.2 - 62.5 %
% Pt, After Reduction With Magnesium	

Specification: PRD.1.ZQ5.10000028590

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.